


Evaluation: I5-XX, 09.11.2021 10:43:51

Curve: I5-XX, 09.11.2021 10:43:51

Sample: 5-XX, 40.8035 mg

Sample Holder: Alumina 70ul
Ceramic

Method: Central South
Forest-Lignin
dt 1.00 s
[1] 30.0..600.0 °C, 30.00 K/min, N2 50.0 ml/min
Sync Enable

Module: TGA/DSC 3+ LF/1100/1024, 12.07.2021 13:32:05

User: METTLER

Curve: I5-XX, 09.11.2021 10:43:51

Sample: 5-XX, 40.8035 mg

Sample Holder: Alumina 70ul
Ceramic

Method: Central South
Forest-Lignin
dt 1.00 s
[1] 30.0..600.0 °C, 30.00 K/min, N2 50.0 ml/min
Sync Enable