

Case number	Slice Thickness	Manufacturer	Model	Contrast enhanced
1	1.0	Toshiba	Aquilion Lightnin	0
2	3.0	GE Medical Syst	Optima CT520 S	0
3	3.0	GE Medical Syst	Revolution Maxir	1
4	5.0	Toshiba	Aquilion	0
5	3.0	GE Medical Syst	Optima CT520 S	0
6	3.0	GE Medical Syst	Revolution EVO	0
7	3.0	GE Medical Syst	Revolution EVO	1
8	5.0	Toshiba	Aquilion Lightnin	0
9	3.0	GE Medical Syst	Revolution EVO	0
10	3.0	GE Medical Syst	Revolution EVO	1
11	3.0	GE Medical Syst	Revolution EVO	0
12	3.0	GE Medical Syst	Revolution Maxir	0
13	3.0	SIEMENS	Emotion 16	0
14	3.0	Toshiba	Aquilion PRIME	0
15	3.0	Siemens Heathir	SOMATOM go.U	0
16	3.0	GE Medical Syst	Optima CT520 S	0
17	3.0	GE Medical Syst	Optima CT520 S	0
18	3.0	GE Medical Syst	Optima CT520 S	0
19	1.0	SIEMENS	Emotion 16	0
20	3.0	GE Medical Syst	Revolution Maxir	0
21	3.0	Toshiba	Aquilion PRIME	1
22	3.0	GE Medical Syst	Optima CT520 S	0
23	5.0	Toshiba	Aquilion Lightnin	1
24	3.0	Toshiba	Aquilion PRIME	1
25	2.0	GE Medical Syst	Revolution Maxir	0
26	3.0	GE Medical Syst	Revolution Maxir	0
27	2.0	GE Medical Syst	Revolution Maxir	0
28	2.0	Toshiba	Aquilion	0
29	3.0	Toshiba	Aquilion PRIME	0
30	3.0	GE Medical Syst	Optima CT520 S	0

31	3.0	Toshiba	Aquilion PRIME	0
32	3.0	GE Medical Syst	Optima CT520 S	0