

Targeting TLR4 with ApTOLL Improves Heart Function in Response to Coronary Ischemia Reperfusion in Pigs Undergoing Acute Myocardial Infarction

Rafael Ramirez-Carracedo ^{1,†}, Laura Tesoro ^{1,†}, Ignacio Hernandez ¹, Javier Diez-Mata ¹, David Piñeiro ², Macarena Hernandez-Jimenez ², Jose Luis Zamorano ³ and Carlos Zaragoza ^{1,*}

Supplemental Material.

Basal cardiac ultrasound parameters of selected pigs included in the study.

PARAMETER	PLACEBO 0	APTOLL 0
IVSD (mm)	10.324±1.672	11.272±1.804
LVIDD (mm)	37.884±5.134	42.237±8.183
LVPWD (mm)	10.629±1.430	11.128±1.537
IVSS (mm)	13.894±1.887	13.991±3.065
LVIDS (mm)	24.371±3.137	26.213±5.726
LVPWS (mm)	15.954±2.648	12.753±1.716
EF (%)	66.79±3.411	68.25±2.571
FS (%)	36.8±3.425	38.2±2.616
HR (bpm)	93.6±12.971	88.8±15.817
CO (ml/sec)	4089.3±343.218	4228.2±470.490