

## **Supplementary Materials**

### *Article*

### **Newly Designed Human-like Collagen to Maximize Sensitive Release of BMP-2 for Markedly Repairing of Bone Defect**

**Zhuoyue Chen<sup>1,2</sup>, Zhen Zhang<sup>1</sup>, Xiaoxuan Ma<sup>2</sup>, Zhiguang Duan<sup>2</sup>, Junfeng Hui<sup>2</sup>, Chenhui Zhu<sup>2</sup>,  
Donggang Zhang<sup>3</sup>, Daidi Fan<sup>2\*</sup>, Lijun Shang<sup>4\*</sup>, Fulin Chen<sup>1\*</sup>**

<sup>1</sup>Provincial Key Laboratory of Biotechnology of Shaanxi; Key Laboratory of Resource Biology and Modern Biotechnology in Western China, Ministry of Education; Northwest University, 229 TaiBai North Road, Xi'an, Shaanxi Province 710069, P.R. China;

<sup>2</sup>Shaanxi Key Laboratory of Degradable Biomedical Materials; Shaanxi R&D Center of Biomaterial and Fermentation Engineering; School of Chemical Engineering, Northwest University, 229 TaiBai North Road, Xi'an, Shaanxi Province 710069, P.R. China;

<sup>3</sup>Yantai Zhenghai Biotechnology Co, Ltd., 10 Hengshan Road, Yantai Economic & Technological Development Area, Yantai, Shandong Province 264000, P.R. China;

<sup>4</sup>School of Chemistry and Biosciences, Faculty of Life Sciences, University of Bradford, Bradford, BD7 1DP, UK;

**\*Correspondence:** [fandaidi@nwu.edu.cn](mailto:fandaidi@nwu.edu.cn); +86-29-88302632 (D. D. Fan); [l.shang1@bradford.ac.uk](mailto:l.shang1@bradford.ac.uk);  
+44-1274-234676 (L. Shang); [chenfl@nwu.edu.cn](mailto:chenfl@nwu.edu.cn); 86029-88302411 (F. L. Chen)

**Figure S1. Mercury porosimetry report of the HLC scaffold** The porosity of the HLC scaffold was over 99%

**Figure S2. Micro-CT images of the rat cranial defect repair** Micro-CT images of control, HLC and HLC-BMP implants groups

**Table S1.** Expression of the recombinant human-like collagens

**Table S2.** Original data of rhBMP-2 Release Kinetics.

**Table S3.** Original data of quantitative estimation of OPN and Runx-2 expression in each implant 2 weeks and 4 weeks after implantation.

**Table S4.** Original data of quantitative estimation of OPN and Runx-2 expression in each implant 8 weeks after implantation.

**Table S5.** Comparison the commercial collagen sponge and the recombinant human-like collagen sponge

**MERCURY POROSIMETRY REPORT**

AutoPore IV 9500 V1.09

Serial: 239

Port: 1/1


Page 9

Sample ID: czy  
Operator:  
Submitter:  
File: C:\9500\DATA\FL-CZY.SMP


LP Analysis Time: 1/5/2017 9:31:39AM  
HP Analysis Time: 1/5/2017 10:15:54AM  
Report Time: 1/5/2017 4:33:12PM

Sample Weight: 0.0292 g  
Correction Type: None  
Show Neg. Int: No

	Cumulative Pore Volume (mL/g)	% of Total Intrusion Volume
	24.1208	99.4510
	24.1208	99.4510
	24.1208	99.4510
	24.1208	99.4510
	24.1208	99.4510
	24.1208	99.4510
	24.1208	99.4510
	24.1208	99.4510
	24.1208	99.4510
	24.1208	99.4510
	24.1208	99.4510
	24.1208	99.4510
	24.1208	99.4510
	24.1208	99.4510
	24.1208	99.4510
	24.1208	99.4510
	24.1208	99.4510
	24.1208	99.4510
	24.1208	99.4510
	24.1208	99.4510
	24.1208	99.4510
	24.1208	99.4510
	24.0192	99.0320



**Figure S1. Mercury porosimetry report of HLC scaffold.** The porosity of the HLC scaffold was over 99%.



**Figure S2. Micro-CT images of the rat cranial defect repair** Micro-CT images of control, HLC and HLC-BMP implants groups

**Table S1. Expression of the recombinant human-like collagens**

Host	Yield	Molecular Wt.	Ref.
<i>E. coli</i> origami (DE3)	260 mg/L in a 10 L fermentor	35 kDa	Tang <i>et al.</i> <sup>1</sup>
<i>E. coli</i> Rosetta (DE3)	Shake flask culture	40 kDa	Yang <i>et al.</i> <sup>2</sup>
<i>E. coli</i> BL21 (DE3)	14 g/L in a 500 L fermentor	93 kDa	Fan <i>et al.</i> <sup>3, 4</sup>

- (1). Tang, Y.; Yang, X.; Hang, B.; Li, J.; Huang, L.; Huang, F.; Xu, Z. Efficient Production of Hydroxylated Human-Like Collagen Via the Co-Expression of Three Key Genes in Escherichia coli Origami (DE3). *Applied Biochemistry & Biotechnology* 2016, 178, 1458-1470.
- (2). Yang, J.; Jie-Ying, Y. U.; Wang, M.; Shao-Min, W. U.; Yi, M. A.; Wang, J. F. Expression, Purification, and Functional Characterization of Recombinant Human-Like Type i Collagen Peptide in E.coli. *Modern Food Science & Technology* 2016.
- (3). Fan, D. D.; Luo, Y.; Mi, Y.; Ma, X. X.; Shang, L. Characteristics of Fed-Batch Cultures of Recombinant Escherichia Coli Containing Human-Like Collagen cDNA at Different Specific Growth Rates. *Biotechnology Letters* 2005, 27, 865-870.
- (4). Zhu, C.; Fan, D.; Duan, Z.; Xue, W.; Shang, L.; Chen, F.; Luo, Y. Initial Investigation of Novel Human-Like Collagen/Chitosan Scaffold for Vascular Tissue Engineering. *Journal of Biomedical Materials Research Part a* 2010, 89A, 829-840.

**Table S2. Original data of rhBMP-2 Release Kinetics.** Day 1~3: Burst release period; Day 3~28 and later: Sustained release period. The percentage of rhBMP-2 release is no difference among loading with different dose rhBMP-2 groups at each time point ( $p > 0.05$ ,  $n=4$ ).

Time (Days)	HLC-BMP 1μg (%)				Average	STD
1	11.75	9.47	7.96	9.78	9.74	1.56
2	13.60	8.66	14.75	14.47	12.87	2.85
3	16.80	19.57	15.58	20.59	18.13	2.34
5	44.77	26.78	34.68	29.88	34.03	7.86
7	35.76	36.80	37.60	42.69	38.21	3.08
14	43.57	40.70	40.68	46.59	42.88	2.82
21	32.85	44.46	56.54	41.73	43.89	9.78
28	40.79	45.69	58.69	33.99	44.79	10.43
Time (Days)	HLC-BMP 2μg (%)				Average	STD
1	15.65	8.57	14.48	9.64	12.09	3.50
2	19.69	15.77	19.58	16.76	17.95	1.99
3	17.54	15.58	23.77	25.80	20.67	4.89
5	46.67	29.68	31.70	39.81	36.96	7.81
7	51.57	39.55	40.76	32.68	41.14	7.81
14	41.80	46.64	42.55	50.57	45.39	4.05
21	45.66	50.46	42.67	44.65	45.86	3.31
28	50.43	43.90	49.12	40.32	45.94	4.69
Time (Days)	HLC-BMP 3μg (%)				Average	STD
1	9.56	8.65	13.77	15.86	11.96	3.43
2	17.85	15.80	10.57	9.68	13.47	3.97
3	23.90	22.75	21.85	17.70	21.55	2.70
5	34.55	36.80	35.46	40.70	36.88	2.71
7	45.60	38.59	34.79	39.90	39.72	4.48
14	45.89	36.78	49.59	42.77	43.76	5.42
21	45.85	40.78	54.86	42.09	45.90	6.35
28	43.87	50.45	44.86	47.98	46.79	3.00
Time (Days)	HLC-BMP 4μg (%)				Average	STD
1	16.58	18.59	9.48	15.47	15.03	3.92
2	16.69	17.69	22.46	17.58	18.61	2.61

3	27.54	18.40	26.50	19.70	23.04	4.65
5	39.35	36.65	39.33	38.54	38.47	1.27
7	42.37	38.16	38.50	39.59	39.66	1.91
14	45.48	49.84	48.27	48.37	47.99	1.82
21	47.80	47.80	49.77	47.68	48.26	1.01
28	51.69	49.78	50.59	39.70	47.94	5.55
Time (Days)	HLC-BMP 5µg (%)				Average	STD
1	16.58	18.59	9.48	15.47	15.03	3.92
2	16.69	17.69	22.46	17.58	18.61	2.61
3	27.54	18.40	26.50	19.70	23.04	4.65
5	39.35	36.65	39.33	38.54	38.47	1.27
7	42.37	38.16	38.50	39.59	39.66	1.91
14	45.48	49.84	48.27	48.37	47.99	1.82
21	45.80	55.80	49.77	43.68	48.76	5.33
28	59.69	39.78	43.59	49.70	48.19	8.69

**Table S3. Original data of quantitative estimation of OPN and Runx-2 expression in each implant 2 weeks and 4 weeks after implantation.**

Time (weeks)/group		OPN (%)			Average	STD
2	Control	0.00	0.00	0.00	0.00	0.00
	HLC	0.00	0.00	0.00	0.00	0.00
	HLC-BMP	21.07	14.13	18.77	17.99	3.54
4	Control	0.00	0.00	0.00	0.00	0.00
	HLC	0.26	0.18	0.24	0.23	0.04
	HLC-BMP	24.77	16.61	22.07	21.15	4.16
Time (weeks)/group		Runx-2 (%)			Average	STD
2	Control	0.00	0.00	0.00	0.00	0.00
	HLC	0.00	0.00	0.00	0.00	0.00
	HLC-BMP	4.97	3.33	4.43	4.25	0.83
4	Control	0.00	0.00	0.00	0.00	0.00
	HLC	0.31	0.21	0.27	0.26	0.05
	HLC-BMP	8.60	5.77	7.67	7.35	1.44

**Table S4. Original data of quantitative estimation of OPN and Runx-2 expression in each implant 8 weeks after implantation.**

Group	OPN (%)			Average	STD
Control	0.09	0.06	0.08	0.08	0.02
HLC- 1µg BMP-2 dimer	35.60	23.87	31.72	30.39	5.98
HLC- 5µg BMP-2 dimer	43.98	29.48	39.18	37.55	7.38
Group	Runx-2 (%)			Average	STD
Control	0.43	0.29	0.38	0.36	0.07
HLC- 1µg BMP-2 dimer	18.66	12.51	16.62	15.93	3.13
HLC- 5µg BMP-2 dimer	23.75	15.92	21.16	20.27	3.99

**Table S5. Comparison the commercial collagen sponge and the recombinant human-like collagen sponge**

Source	Company	Release rate	Effect
bovine achilles tendon -extracted collagen <sup>5</sup>	Helistat®, Integra Life Sciences Corporation , USA	5 µg BMP-2 dose group released $55.46 \pm 2.01\%$ within 1 hour, $84.94 \pm 8.68\%$ within 1 day, and 100% at day 2. The 1 µg BMP-2 group released $34.73 \pm 8.25\%$ within 6 hours, and 45% at the third day.	The collagen sponge containing 1 µg BMP-2 formed the new bone volume was $0.6 \pm 0.3 \text{ mm}^3$ at 8 weeks post-implantation into the rat muscle bag.
recombinant from <i>E. coli</i>	Our lab.	HLC containing 5 µg BMP-2 gel released $15.03 \pm 3.92\%$ of BMP-2 at day 1, $18.61 \pm 2.61\%$ at day 2 and $48.19 \pm 8.69\%$ at day 28. The HLC containing 1 µg BMP-2 released $9.74 \pm 1.56\%$ at day 1, $18.13 \pm 2.34\%$ at day 3, and $44.79 \pm 10.43\%$ at day 28.	Using the similar rat cranial defect repair model, the repair rate is $65.21 \pm 4.44\%$ at 2 weeks postoperatively. At 4 weeks postoperatively, cranial defect repair rate of the designed HLC loaded with 1 µg BMP-2 group achieved $79.95 \pm 7.57\%$ .
bovine skin -extracted collagen <sup>6</sup>	I-PC, Koken Co. Ltd. Tokyo, Japan	Collagen gel containing 60 µg of BMP-2 released 6% of BMP-2 at day 1. About 33% of the total BMP-2 was released from the collagen gel within 15 days.	None of the new bone tissue appeared in the mouse cranial defect area at 2 weeks postoperatively. At 4 weeks postoperatively, the cranial defect repair rate of the commercial sponge loaded with almost 1 µg BMP-2 group was $40.24 \pm 9.70\%$ .

(5). Bhakta, G.; Lim, Z. X. H.; Rai, B.; Lin, T.; Hui, J. H.; Prestwich, G. D.; Wijnen, A. J. V.; Nurcombe, V.; Cool, S. M. The Influence of Collagen and Hyaluronan Matrices On the Delivery and Bioactivity of Bone Morphogenetic Protein-2 and Ectopic Bone Formation. *Acta Biomaterialia* 2013, 9, 9098-9106.

(6). Zhao, J.; Shinkai, M.; Takezawa, T.; Ohba, S.; Chung, U.; Nagamune, T. Bone Regeneration Using Collagen Type I Vitrigel with Bone Morphogenetic Protein-2. *Journal of Bioscience and Bioengineering* 2009, 107, 318-323.