

Article

# Gold Nanoparticles/Carbon Nanotubes and Gold Nanoporous as Novel Electrochemical Platforms for L-Ascorbic Acid Detection: Comparative Performance and Application

Cristina Tortolini <sup>1</sup>, Federico Tasca <sup>2</sup>, Mary Anna Venneri <sup>1</sup>, Cinzia Marchese <sup>1</sup> and Riccarda Antiochia <sup>3,\*</sup>

<sup>1</sup> Department of Experimental Medicine, Sapienza University of Rome, Viale Regina Elena 324, 00166 Rome, Italy; cristina.tortolini@gmail.com (C.T.); maryanna.venneri@uniroma1.it (M.A.V.); cinzia.marchese@uniroma1.it (C.M.)

<sup>2</sup> Departamento de Química de los Materiales, Facultad de Química y Biología, Universidad de Santiago de Chile, Casilla 40, Correo 33, Sucursal Matucana, Santiago 9170022, Chile; federico.tasca@usach.cl

<sup>3</sup> Department of Chemistry and Drug Technologies, Sapienza University of Rome, P.le Aldo Moro 5, 00185 Rome, Italy

\* Correspondence: riccarda.antiochia@uniroma1.it

**Citation:** Tortolini, C.; Tasca, F.;

Venneri, M.A.; Marchese, C.;

Antiochia, R. Gold

Nanoparticles/Carbon Nanotubes

and Gold Nanoporous as Novel

Electrochemical Platforms for L-

Ascorbic Acid Detection: Comparative

Performance and Application.

*Chemosensors* **2021**, *9*, 229.

<https://doi.org/10.3390/chemosensors9080229>

Academic Editor: Dale A. C.

Brownson and Alexey Glushenkov

Received: 26 June 2021

Accepted: 3 August 2021

Published: 16 August 2021

**Publisher's Note:** MDPI stays neutral with regard to jurisdictional claims in published maps and institutional affiliations.



**Copyright:** © 2021 by the authors.

Licensee MDPI, Basel, Switzerland.

This article is an open access article


distributed under the terms and

conditions of the Creative Commons

Attribution (CC BY) license

([http://creativecommons.org/licenses](http://creativecommons.org/licenses/by/4.0/)

[/by/4.0/](http://creativecommons.org/licenses/by/4.0/)).



**Figure S1.** DP voltammograms recorded in 0.1 M PBS pH = 7.4, KCl = 0.1 M at different concentrations of AA (10, 150, 250, 400, 800  $\mu\text{M}$ ) using h-nPG electrode. Inset: calibration curve of AA.