

Table S1. Liquid media composition.

Medium	Composition (per Liter basis)
Cell Suspension/ Cellulolytic¹	<ul style="list-style-type: none"> • 240 mg KH₂PO₄ • 240 mg K₂HPO₄ • 480 mg (NH₄)₂SO₄ • 480 mg NaCl • 64 mg CaCl₂ • 2H₂O • 100 mg MgSO₄ • 7H₂O • 600 mg cysteine hydrochloride • 1,000 mg Trypticase • 500 mg yeast extract • 1 mg phenylacetic acid • 3.1 mL VFA solution (7.15 mol L-1 acetic acid, 1.92 mol L-1, propionic acid, 783.66 mmol l-1 butyric acid, 222.28 mmol l-1 valeric acid, 216.19 mmol l-1 isovaleric acid, 264.13 mmol l-1 isobutyric acid, 227.87 mmol l-1 2-methylbutyric acid).
Hyper-Ammonia Producing Bacteria²	<ul style="list-style-type: none"> • 240 mg KH₂PO₄ • 240 mg K₂HPO₄ • 480 mg Na₂SO₄ • 480 mg NaCl • 64 mg CaCl₂ • 2H₂O • 100 mg MgSO₄ • 7H₂O • 600 mg cysteine hydrochloride • 10 mL vitamin solution (4.14.75 μmol L-1 pyridoxamine 2 HC., 531.41 μmol L-1 riboflavin, 592.99 μmol L-1, 592.99 μmol L-1 thiamine HCl, 1.64 mmol l-1 nicotinamide, 419.73 μmol L-1 CaD pantothenate, 484.66 μmol L-1 lipoic acid, 72.92 μmol L-1 p-aminobenzoic acid, 11.33 μmol-1 folic acid, 20.47 μmol L-1 biotin, 3.69 μmol L-1 cobalamine, 491.11 μmol L-1 pyridoxal HCl, 591.09 μmol L-1 pyridoxine) • 5 mL trace mineral solution (5.44 mmol l-1 Na₄ EDTA, 719 μmol L-1 FeSO₄, 34.78 μmol ZnSO₄, 1.01 mmol l-1 MnCl₂, 323.47 μmol L-1, H₃BO₃, 96.61 μmol L-1 CoCl₂. 5.87 μmol L-1 CuCl₂, 8.41 μmol L-1 NiCl₂ 6.98 μmol L-1 Na₂MoO₄)
Lactate-Utilizing Bacteria³	<ul style="list-style-type: none"> • 352 mg KH₂PO₄ • 352 mg K₂HPO₄ • 704 mg Na₂SO₄ • 701 mg NaCl • 94 mg CaCl₂ • 2H₂O • 147 mg MgSO₄ • 7H₂O • 100 mg cysteine • 220 mmol lactic acid • 7.67 μmol hemin • 10 mL vitamin solution (as described above) • 10 mL mineral (as described above) • 10 mL VFA (as described above)

¹Prepared based on protocol described by Stack et al. [16]. ²Prepared as previously described by Russell et al. [18] and Chen and Russell [19]. ³Prepared as described by Mackie and Heath [20].