Anton Paar GmbH
Anton-Paar-Straße 20
8054 Graz
Austria

Anton Paar Monowave

Serial Number: 81247851
Instrument Software Version: 3.20.7338.3

Processing Protocol


- Experiment Name: sanja_copy13
- Experiment Date: 9/6/2023 2:33:35 PM
- User: Dule

Steps

Step	Program	Temperature	Time	Power	Cooling	Stirrer Speed
		°C	hh:mm:ss	W		rpm
1	Constant power	119	-	50	On	600
2	Constant power	120	-	50	On	600
3	Cool down	55	-	-	On	600

Experiment Result

- Result: OK
- User name of last IR sensor adjustment: Administrator
- Vial Type: Glass vial G30



Raw data is accessible for download from the "Browse Results" dialog on Monowave 300 (see serial number above)

Wednesday, September 06, 2023 _____
(signature)

Figure S4. Example of a power-controlled mode at 50W for 2 min, with a graphical representation of the temperature change during the experiment