Supplemental Data 3: Negative controls for Figure 3 (in main text)

Immunofluorescence staining of first trimester placental villi (first row) and decidua basalis (second row) from placental explant cultures (9th week of pregnancy) using Isotype control antibodies. Left panel: DAPI staining to visualize nuclei. Middle panel: Isotype Control for DUSP9 and PHLPP with anti-rabbit 488 secondary antibody (green). Right panel: Isotype control for HLA-G in combination with anti-mouse 568 secondary antibody (red). Bars indicate 50 µm. DB, decidua basalis; pV, placental villus;

