Table S2. In vitro medium formulations

Reagent	Volume/Mass
Sf-900™ II SFM (Thermo Fisher Scientific, New York, NY, USA)	10 mL
Penicillin-Streptomycin Solution，100× (Servicebio, Wuhan, China)	100 μL
Glutathione, Reduced (Shanghai RYON Biotechnology Co., Ltd,Shanghai, China)	9 mg
Kanamycin	20 μL
Ampicillin	20 μL
JH-Ⅲ/20E (Shanghai Yuanye Bio-Technology Co., Ltd., Shanghai, China)	final concentration to 5 μM